




Technowood Radial Access Essentials



Hydrophobic surface and wet-cloth cleanable.
 Cleanable with alcohol-based solution.
 Arm Rest (urethane-based) is elastic and absorbs impact.
 Adjustable holding grip (urethane-based)
 Arm Rest attaches to Platform conveniently with velcro tape.

Simply insert the platform **A** between operating table and mattress for an easy setup.
 Platform (PTE-based) provides excellent transparency.
 Platform is slightly angled for patient's comfort as well as easy puncture site management.



<u>Product</u>	<u>Description</u>	<u>Height</u>	<u>Length</u>	<u>Width</u>	<u>Remarks</u>
Arm Rest	left & right	19 cm (7.48")	35 cm (13.78")	12 cm (4.72")	holding grip (adjustable)
Platform	left & right				

<u>Cat. #</u>	<u>Specifications</u>
EX-FT-TEOKI-R	Arm Rest (right) + Platform (right)
EX-FT-TEOKI-L	Arm Rest (left) + Platform (left)

★ **JAPAN**
Technowood International Corporation
 30-19 Kohoku, 4-Chome
 Adachi-Ku, Tokyo 123-0872
 Japan

★ **USA**
Technowood America Corporation
 3001 Red Hill Avenue
 Building 2, Suite 201
 Costa Mesa, CA 92626
 U.S.A.

★ **Philippines**
Technowood Manufacturing Corporation
 N2840 Bayanihan Street Extension
 CSEZ Clarkfield, Pampanga 2009
 Philippines

Tel **81-3-3898-5252**
 Fax: **81-3-3898-5252**
 Web-Site <http://www.technowood.co.jp>
 Contact info-tic@technowood.co.jp

Tel **1-714-434-8713**
 Fax: **1-714-434-8715**
 Web-Site <http://www.technowood.com>
 Contact info@technowood.com

Tel **63-45-599-5570 or 63-45-599-5571**
 Fax: **63-45-599-5574**
 Web-Site <http://www.technowood.com.ph>
 Contact info@technowood.com.ph