

Model No. 801-100



### Innovation in cardiac catheterization

- Compact Design
- Simple Operation
- Precise Injection
- Reduce Fatigue

### Performance

- Injection power generated by pedestal operation.
- Significant reduction in fatigue over the current angiography by hand operation.
- Ideal for small size catheters.
- Reduced contrast media volume.
- No need for large volume power injector syringe.
- Accommodate 10 cc & 20 cc control syringe.

### Mobility

- All in one system can be easily moved between rooms.

### Space Saving

- Very small base size never interferes your procedural efficacy in the cath lab.

### Simplicity

- Easy operation.
- Minimum maintenance - save maintenance cost.


### Comprehensive

- Increased efficiency and productivity and patient management.
- Suitable for Coronary Angiography, LV Angiography, Aorta Flush.





with sterile cover (cat. # 801-101)


Technowood Foot Injector is designed to optimize the cardiac catheterization with simplicity, efficiency, and productivity while minimizing the cost associated with cath lab procedures.



10 cc control syringe



20 cc control syringe



\* Above test data at our lab provides the vital information that the Technowood Foot Injector can accomplish the coronary angiography with a 10 cc control syringe and LV angiography, aorta flush with a 20 cc control syringe compared with the conventional 30 cc or greater syringes for the power injector, which further helps the cath lab save cost on the contrast media.

Description	Model #	Specifications	Feature	Accessories
Foot Injector	801-100	W: 14" L: 28" H: 50" 33 lbs	Movable Syringe Arm Holder Adjustable Mast Height Swiveled Arm Syringe Flow Rate Setting Contrast Media Remaining Volume Setting	Disposable Control Syringe (10 ml & 20 ml) Disposable Sterile Cover (Cat. # 801-101)

- Specifications are subject to change without notice.
- Technowood Foot Injector

Japanese M.H.W. Approval # 21300BZZ00613000

Manufactured By:  
**Tonokura Ika Kogyo Co., Ltd.**  
 1-13 Hongo, 5-Chome  
 Bunkyo-Ku, Tokyo, 113-8435  
 Japan

TEL 03-3813-7751  
 FAX 03-3813-8537

U.S. Office:  
**Technowood America Corporation**  
 1360 Reynolds Ave.  
 Suite 112  
 Irvine, CA 92614  
 U.S.A.

TEL 949-794-0640  
 FAX 949-794-0642  
 Web Site <http://www.technowood.com>  
 E-mail [info@technowood.com](mailto:info@technowood.com)